

# RAPID FAST FIELD-CYCLING MRI USING KEYHOLE IMAGING

Peter James Ross<sup>1</sup> and David J. Lurie<sup>1</sup>

<sup>1</sup>Aberdeen Biomedical Imaging Centre, University of Aberdeen, Aberdeen, Aberdeen City, United Kingdom

## Purpose

Fast Field-Cycling MRI (FFC-MRI)<sup>1</sup> is an emerging technique that adds a new dimension to conventional MRI by making it possible to rapidly vary  $B_0$  during a pulse sequence. By doing this it is possible to observe how the NMR relaxation rates of biological tissues vary with magnetic field strength - information which can be employed as a useful contrast mechanism. To date we have used FFC-MRI to perform spatially-selective relaxometry using an adapted PRESS sequence<sup>2</sup> and also relaxometric imaging (a set of  $R_1$  images at a range of field strengths) using a gradient echo sequence. Relaxometric imaging collects more information than selected-volume relaxometry but its application is limited by lengthy scan times, since the entirety of k-space is acquired at each field strength. For high-resolution imaging, or where images are collected at multiple evolution fields, scan times can become unacceptably long.

In this work we have made use of the keyhole MRI technique<sup>3</sup> in order to speed up FFC-MRI. By collecting data for the whole of k-space at the beginning of each scan and thereafter only updating the low spatial-frequency region of k-space with each subsequent field-cycling experiment, contrast derived from the FFC technique is maintained while the scan time is dramatically reduced.

## Methods

Imaging was carried out on a home-built, whole-body, field-cycling imager with a 59 mT detection field and a coaxial resistive offset magnet which provides field-cycling capabilities<sup>4</sup>. The system uses a commercial console (MR Solutions, U.K.).

For each experiment, reference saturation-recovery and inversion-recovery images were acquired at the detection field using a conventional spin-echo sequence. Following this, for each magnetic field of interest, the "keyhole" portion of k-space - corresponding to the central 25% of k-space was acquired using a field-cycling inversion-recovery spin-echo sequence. During the inversion recovery period  $B_0$  was rapidly switched to a different field, referred to as the evolution field, and  $M_z$  was allowed to relax at that field for an evolution period, typically of the order of  $T_1$ . The field was then switched back to the detection field and the imaging sequence was performed.

The data were then reconstructed by combining the initially-acquired high spatial-frequency part of k-space with the keyhole central portion of k-space, obtained at each evolution field value. In this way a full-resolution image was constructed at each evolution field value, requiring a scan time of 25% compared to conventional imaging. The technique was used to derive  $R_1$  dispersion curves ( $R_1$  vs Larmor frequency) from a phantom consisting of cross-linked bovine serum albumin (BSA) across a field range of 1.4 MHz to 2.5 MHz (32 mT to 59 mT).

Brain images were also collected from a volunteer (Fig 2A) in order to validate that the method generated artefact free images *in vivo*. Images were collected for evolution fields of 49 mT and 59 mT (proton Larmor frequencies of 2.1 MHz and 2.5 MHz) and from these images  $R_1$  maps were derived. All data processing and image reconstruction was performed using in-house software developed using MATLAB R2014a.

## Results and Discussion

$R_1$  dispersion curves derived from BSA using the keyhole technique show excellent agreement with results obtained using a conventional full k-space scan. (Fig.1). The results show a distinct increase in  $R_1$  at close to 2.1 MHz (49 mT). This feature is known as a quadrupole peak, and arises due to  $^1\text{H}$ - $^{14}\text{N}$  cross relaxation effects, which occur at well-known field strengths. The high  $^{14}\text{N}$  concentration typically present in protein rich tissue which leads to these peaks can be exploited to generate protein sensitive contrast in FFC-MRI images.

The  $R_1$  maps (Fig. 2B and 2C) generated from brain images show an increase in  $R_1$  at 49 mT, which was used to generate a  $\Delta R_1$  map (Fig. 2D) which shows a higher increase in  $R_1$  in predominately white matter regions.

## Conclusions

This work has demonstrated that the keyhole technique can readily be applied to FFC-MRI and used to obtain a 4-fold or greater speed up in scan times while still retaining the same contrast as standard FFC-MRI methods. The rich  $\Delta R_1$  contrast present in brain images indicates that FFC-MRI has potential application in the characterization of neurodegenerative conditions where subtle changes in  $R_1$ , which may not be visible on conventional  $T_1$ -weighted imaging could be used as an early marker of disease.

The reduction in scan time achieved by use of the keyhole technique will significantly improve the applicability of FFC-MRI in volunteer and clinical studies, which we are currently working towards.

## Acknowledgements

The author acknowledges funding from the EPSRC through the Centre for Doctoral Training in Integrated Magnetic Resonance.

## References

- [1] Lurie D.J. *et al.* Fast field-cycling magnetic resonance imaging. *Comptes Rendus Phys.* 2010;11:136-148.
- [2] Pine K.J., Davies G.R., Lurie D.J. Field-cycling NMR relaxometry with spatial selection. *Magn. Reson. Med.* 2010;63:1698-702
- [3] van Vaals J.J. *et al.* "Keyhole" method for accelerating imaging of contrast agent uptake. *J. Magn. Reson. Imaging* 3:671-5.
- [4] Lurie D.J., Foster M.A., Yeung D., Hutchison J.M. Design, construction and use of a large-sample field-cycled PEDRI imager. *Phys. Med. Biol.* 1998;43:1877-86.




Figure 2. A: Brain image collected from a volunteer at 59 mT  
B: Keyhole  $R_1$  map collected at 59 mT  
C: Keyhole  $R_1$  map collected at 49 mT  
D:  $\Delta R_1$  map generated from the subtraction of B from C.  
Image parameters were: matrix size = 128x128, FOV = 280 mm, THK = 15 mm, NEX = 6, TE = 30 ms, TR = 1500 ms, evolution time = 200 ms, field-cycling ramp time = 30 ms.




Figure 1: Dispersion curves for a phantom of cross-linked BSA obtained using the keyhole method (red diamonds) show good agreement with results obtained using a conventional field-cycling spin-echo sequence (black circles).